purpose.

MR SAFETY INFORMATION

The MetaFix™ Hip stem has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the MetaFix™ Hip Stem in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

ADVERSE REACTIONS

Complications can potentially occur following any joint replacement surgery. The following complications have occurred in some patients after their hip joint replacement surgery.

- Hip ball (femoral head) and hip socket (acetabulum) may separate (hip dislocation);
- Device loosening from the surrounding bone;
- Allergic reaction to the implant's materials;
- Audible sounds during motion;
- Premature wear or breakage of the implants;
- Bone loss around the implant;
- Periprosthetic fracture;
- Change in the length of the treated leg;
- Hip pain and stiffness;
- Loss of hip flexibility of the hip joint;
- Nerve damage;
- Calcification.

STERILITY

The MetaFix™ Hip stem implants are supplied sterile. Check the integrity of the packaging

carefully. Do not use if the packaging is open or damaged. Do not sterilize. For single use only. The product is not labelled "pyrogen free".

INSTRUMENTS

Instruments used to implant the device are supplied non-sterile and must be sterilised prior to use. Instruments should be cleaned according to the Corin cleaning instructions supplied with the instrument set. After cleaning and prior to sterilisation, the reusable device(s) should be double-wrapped or packaged in FDA cleared (USA only) wraps or pouches. Wrapping should be performed using the appropriate wrapping method (e.g. AAMI CSR wrapping technique). The following sterilisation method has been validated, based on AAMI/ANSI/ISO guidelines and recommendations:

- Method: Moist-Heat Sterilization
- Cycle: Pre-Vacuum (Pre-Vac)
- Temperature: 270 F (132 C)
- Exposure Time: 4 minutes
- Pressure: 2-15 PSIA
- Dry-Time: 30 minutes (minimum, in chamber)
- Cool-Time: 60 minutes (minimum, at room temperature)

In the USA other sterilization methods for instruments should not be used as they have not been validated by Corin.

After sterilizing, the reusable device(s) should be left on the sterilizer cart, untouched, for a minimum of one hour at room temperature, or until adequately cooled for safe handling. Other sterilization

methods and cycles may also be suitable, however, individuals or hospitals are advised to validate whichever method they deem appropriate at their institution. Ethylene Oxide (EtO) sterilization and cold sterilization techniques are not recommended.

remove surgical debris prior to closure. Fracture, particularly of smaller sized stems, is most likely to occur in patients who are young, physically active, and/or heavy. The following conditions, either singly or in combination, can put the patient at higher risk of failure of their joint replacement due to severe loading of the affected extremity: obesity or excessive patient weight, manual labour, active sports participation, high levels of patient activity, likelihood of falls, alcohol or drug addiction, and other disabilities as applicable.

Surgeons should give consideration to bone and muscle quality. Clean gloves should be worn when handling implants. Avoid scratching or denting implant surfaces. Do not modify implants. Do not reuse implants, an implant may appear undamaged but previous stress may have created imperfections that would reduce the service life of the implant. Surgeons should be thoroughly familiar with the Corin MetaFix™ Hip surgical technique, the implants, and the instruments prior to performing the surgery.

Corin designs specialized instruments for their joint replacement systems to aid in the accurate implantation of the prostheses. The use of instruments or implant components from other systems can result in inaccurate fit, sizing, and device failure. Intra-operative fracture or breaking of instruments has been reported. Surgical instruments are subject to wear with normal usage. Instruments which have experienced extensive use or excessive force, are susceptible to fracture. Surgical instruments should only be used for their intended

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