



OMNIlife science, Inc.

480 Paramount Dr.
Raynham, MA 02767
www.omnils.com
(508)824-2444

The OMNI Hip™ System K1 Hip Stem**

SYMBOLS Glossary per ISO 15223-1

	Medical Device Manufacturer
	Date of Manufacture
	Use By Date
	Do not Re-use
	See Instructions for Use
	Do Not Use if Package is Damaged
	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.
QTY	Quantity
	Sterilized Using Ethylene Oxide
	Batch Code
	Catalogue Number
	Medical Device
	MRI Safety Information

PRODUCT HANDLING

Implants are provided sterile and should always be stored unopened in their respective protective containers. Prior to use, inspect package for damage which may compromise sterility. If packaging has been opened or damaged, contact manufacturer's representative. When unpacking the implant, verify the labeling for correct Cat. No. and size. When removing the implant from its packaging, the sterile technique must be observed. Protect the implant from contact with objects that may damage the surface finish. Inspect each implant prior to use for visual damage. Procedures

for implanting and removal are available upon request.

DESCRIPTION

The OMNI K1™ Hip stem is a monoblock press-fit hip stem for uncemented applications. It is available in sizes 0 to 12 with various neck options. The stem surface is grit blasted with a proximal CP titanium plasma spray. The OMNI K1 Hip Stem is for use with OMNI and Corin modular femoral heads.

- Use of heads with a shorter (more negative) offset than -4mm with any K1 stem is not indicated.

For more details regarding the compatibility between implants, please contact your OMNI representative or you could find more details on the following link www.coringroup.com/compatibility.

MATERIAL

- Titanium alloy (ASTM F 136) and unalloyed titanium (ASTM F 1580).

INDICATIONS FOR USE

The indications for use of the K1 Hip Stem in hip arthroplasty include the following conditions, as appropriate:

- Non-inflammatory degenerative joint disease, including osteoarthritis and avascular necrosis;
- Rheumatoid arthritis;
- Correction of functional deformity;
- Congenital dislocation;
- Revision procedures where other treatments or devices have failed;
- Femoral neck and trochanteric fractures of the proximal femur.

The K1 Hip Stem is indicated for cementless use only and single use implantation.

CONTRAINDICATIONS

Absolute contraindications include:

- Infection or sepsis or osteomyelitis;
- Insufficient bone structure or quality which may affect the stability of the implant;
- Rapid joint destruction or bone absorption;
- Skeletal immaturity;
- Muscular, ligamentous, neurological, vascular deficiencies or poor skin coverage, which may compromise the affected extremity;

- Alcoholism or other addictions;
- Sensitivity to the implant materials;
- High levels of physical activity (e.g. competitive sports, heavy physical labor);
- Obesity that can produce loads on the prosthesis, which can lead to fixation failure or prosthesis breakage or fracture.

Relative contraindications include:

- Uncooperative patient or a patient with neurological disorders and incapable of following instruction;
- Metabolic disorders which may impair bone formation or bone quality;
- Distant foci of infections.

WARNINGS AND PRECAUTIONS

While these implants are not intended to withstand activity levels and loads of normal healthy bone, they are a means of restoring mobility and reducing pain for many patients.

In using total joint implants, the surgeon should be aware of the following:

- The correct selection of implant components is extremely important. The potential for success in total joint replacement is increased by the selection of the proper size, shape and design of the implant. Total joint prostheses require careful seating and adequate bone support, and should be restricted to limited functional stress. The surgeon is to be thoroughly familiar with the implant, instruments, and surgical procedure prior to performing surgery.
- In selecting patients for joint replacement surgery, the following factors can be of extreme importance to the eventual success of the procedure:
 1. The patient's weight. An overweight or obese patient can produce loads on the prosthesis, which can lead to failure of the prosthesis. This becomes a major consideration when a small prosthesis must be used.
 2. The patient's occupation or activity. If the patient is involved in an occupation or activity, which includes substantial walking, running, lifting, or muscle strain, the resultant forces can cause failure of the fixation, the device or both.
 3. A condition of senility, mental illness or alcoholism. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions, leading to failure or other complications.
 4. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be performed prior to material selection or implantation.
 5. Certain degenerative diseases. In some cases, the progression of degenerative disease may be so

advanced at the time of implantation that it may substantially decrease the expected useful life of the appliance. For such cases, total joint replacement can only be considered a delaying technique or temporary relief.

- The correct handling of the implant is extremely important. Care must be taken to protect surfaces from nicks and scratches which could become the focal point for failure. Do not tamper with the implant as contouring or bending of the implant may reduce its service life and may cause immediate or eventual failure under load.
- Standard cleaning procedures cannot be relied upon to remove contamination from porous coating.
- A surgical implant should not be reused. Even though a used implant may appear undamaged, it may have small defects and internal stress patterns, which may lead to failure. Use only new prosthesis of the current design.
- Resterilization of the device is not recommended.
- Bone excision should be limited to the amount necessary to accommodate the implants. Prior to closure, the surgical site should be thoroughly cleaned of bone chips, bone cement or other detritus that may cause a third body wear problem. Range of motion should be checked for impingement or instability.
- Postoperative care is important. The patient should be instructed on the limitations of these devices and should be cautioned regarding load-bearing, ranges of motion, and activity levels permissible. Excessive physical activity and trauma affecting the replaced joint have been implicated in premature failure by loosening, fracture and/or wear of the prosthesis implant. Early load-bearing should be carefully controlled. The patient should be advised to report any related pain, decrease in range of motion, swelling, fever, and unusual incidences.
- The modular head must be firmly seated to prevent disassociation. Scratching of modular heads and tapers should be avoided. Repeated assembly and disassembly of the head components could compromise a critical locking action. The head should be changed only when

clinically necessary. The interfaces should be clean and free from debris prior to assembly.

POSSIBLE ADVERSE EFFECTS

The possible adverse effects are similar to those occurring with any hip arthroplasty and include the following:

- Dislocation or subluxation of the hip ball (femoral head) and hip socket (acetabulum).
- Loosening or migration of components due to trauma and/or loss of fixation.
- Accelerated wear of the polyethylene articulating surfaces of acetabular components. Such wear may be initiated by particles of cement, metal, or other debris which can cause abrasion of the articulating surfaces. Accelerated wear shortens the useful life of the prostheses and may lead to early revision surgery to replace the worn components.
- Histiocytic granuloma formation and osteolysis around the implant due to wear debris.
- Aseptic, lymphocyte-dominated vasculitis-associated lesion (or pseudotumor) formation due to fretting wear and/or corrosion at modular implant junction(s).
- Fracture of the implant as the result of strenuous activity, obesity, improper alignment, inadequate fixation or extreme duration of service.
- Urological complications, especially urinary retention and infection.
- Dislocation, wear, dissociation, or fracture due to device impingement during hip articulation.
- Other complications associated with general surgery, such as drugs, ancillary devices, and blood loss or transfusions.

Intraoperative and early postoperative complications can include:

- Damage to blood vessels;
- Temporary or permanent neuropathies;
- Traumatic arthrosis of the knee from Intraoperative positioning of the extremity;
- Cardiovascular disorders including venous thrombosis, pulmonary embolism, or myocardial infarction;
- Hematoma;

- Delayed wound healing;
- Infection;
- Femoral perforation;
- Fracture of the femur;
- Undesirable shortening or lengthening of the limb.

Late postoperative complications can include:

- Aggravated problems of the knee or ankle of the affected limb or contralateral extremity by leg length discrepancy, too much femoral medialization, or muscle deficiency;
- Femoral fracture by trauma or excessive loading, particularly in the presence of poor bone stock;
- Periarticular calcification or ossification, with or without impediment to joint mobility;
- Inadequate range of motion due to improper selection or positioning of components, by femoral impingement and periarticular calcification;
- Excessive joint pressures and pain with ambulation due to excessive scarring of the joint capsule and surrounding tissues;
- Infection;
- Trochanteric non-union due to inadequate reattachment and/or early weight bearing;
- Trochanteric avulsion as a result of excessive muscular weakening.

CAUTION

Disposal of implants should be carried out using the hospital's standard method for non-biodegradable non-combustible medical waste.

MRI SAFETY INFORMATION

The implants have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the implants in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

All rights reserved. OMNIHIP™, OMNI Interface™ Acetabular System, OMNI Mod™ and OMNI K1™ are trademarks of OMNIlife science, Inc. Additional information about the OMNI Hip System may be obtained from OMNIlife science, Inc.