GB Attention Operating Surgeon

IMPORTANT Please read carefully before implanting this product

DESCRIPTION

The Corin TriFit TS Hip is a double tapered blade stem manufactured from Ti6Al4V Titanium alloy with a layer of commercially pure Titanium and electrochemically deposited calcium phosphate coating applied. The TriFit TS Hip is available in standard and lateralized offsets in 127° CCD angle. The device is intended to be used with Corin 12/14 modular taper heads.

The TriFit TS Hip is intended to provide increased patient mobility and reduce pain by replacing the damaged hip joint articulation in patients where there is evidence of sufficient sound bone to seat and support the components.

INDICATIONS

The indications for the TriFit TS Hip as a total hip arthroplasty, and when used in combination with a Corin hemi-arthroplasty head, as a hip hemi-arthroplasty, include:

- Non-inflammatory degenerative joint disease including osteoarthritis and avascular necrosis
- Rheumatoid arthritis
- Correcion of functional deformity
- Treatment of non-union and femoral neck fractures
- Developmental Dysplasia of the Hip (DDH) and

Corin TriFit TS Hip, showing the modular taper head and stem.

 Congenital Dysplasia of the Hip (CDH).

The TriFit TS Hip is indicated for cementless use only.

CONTRAINDICATIONS

- Active Infection
- Osteomyelitis
- Poor bone quality
- Marked bone loss or bone resorption
- Metabolic disorders which may impair bone formation
- Vascular insufficiency
- Muscular atrophy or neuromuscular disease
- Allergy to implant material
- Uncorrectable deformity
- Do not use in combination with components from non Corin implant systems

WARNINGS AND PRECAUTIONS

Failure of the implant may result from: selection of improper components, mal-alignment or mal-positioning of the components, inadequate fixation of components, failure to remove surgical debris prior to closure. Fracture, particularly of smaller sized stems, is most likely to occur in patients who are young, physically active, and/ or heavy. The following conditions, either singly or in combination, can put the patient at higher risk of failure of their joint replacement due to severe loading of the affected extremity: obesity or excessive patient weight, manual labour, active sports participation, high levels of patient activity, likelihood of falls, alcohol or drug addiction, and

Corin TriFit TS Hip, showing the modular taper head and stem.

other disabilities as applicable.

Surgeons should give consideration to bone and muscle quality. Clean gloves should be worn when handling implants. Avoid scratching or denting implant surfaces. Do not modify implants. Do not reuse implants; an implant may appear undamaged but previous stress may have created imperfections that would reduce the service life of the implant. Surgeons should be thoroughly familiar with the TriFit TS Hip surgical technique, the implants, and the instruments prior to performing the surgery.

Corin designs specialised instruments for their joint replacement systems to aid in the accurate implantation of the prostheses. The use of instruments or implant components from other systems can result in inaccurate fit, sizing, and device failure. Intra-operative fracture or breaking of instruments has been reported. Surgical instruments are subject to wear with normal usage. Instruments which have experienced extensive use or excessive force, are susceptible to fracture. Surgical instruments should only be used for their intended purpose.

MR SAFETY INFORMATION

The TriFit TS stem has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the TriFit TS Stem in the MR environment is unknown. Scanning a patient who has this device

may result in patient injury.

ADVERSE REACTIONS

Complications can potentially occur following any joint replacement surgery. The following complications have occurred in some patients after their hip joint replacement surgery.

- Hip ball (femoral head) and hip socket (Acetabulum) may separate (hip dislocation);
- Device loosening from the surrounding bone;
- Allergic reaction to the implant’s materials;
- Audible sounds during motion;
- Premature wear or breakage of the implants;
- Bone loss around the implant;
- Periprosthetic fracture
- Change in the length of the treated leg;
- Hip pain and stiffness;
- Loss of hip flexibility of the hip joint;
- Nerve damage;
- Calcification.

STERILITY

The TriFit TS Hip implants are supplied sterile. Check the integrity of the packaging carefully. Do not use if the packaging is open or damaged. Do not resterilise. For single use only. The product is not labelled “pyrogen free”.

INSTRUMENTS

Instrments used to implant the device are supplied non-sterile and must be sterilised prior to use. Instruments should be cleaned according

Corin TriFit TS Hip, showing the modular taper head and stem.

to the Corin cleaning instructions supplied with the instrument set. After cleaning and prior to sterilisation, the reusable device(s) should be double-wrapped or packaged in FDA cleared (USA only) wraps or pouches. Wrapping should be performed using the appropriate wrapping method (e.g. AAMI CSR wrapping technique). The following sterilisation method has been validated, based on AAMI/ANSI/ISO guidelines and recommendations:

- Method: Moist-Heat Sterilisation
- Cycle: Pre-Vacuum (Pre-Vac)
- Temperature: 270°F (132°C)
- Exposure Time: 4 minutes
- Pressure: 28.5 PSIG
- Dry-Time: 30 minutes (minimum, in chamber)
- Cool-Time: 60 minutes (minimum, at room temperature)

In the USA other sterilization methods for instruments should not be used as they have not been validated by Corin.

After sterilising, the reusable device(s) should be left on the steriliser cart, untouched, for a minimum of one hour at room temperature, or until adequately cooled for safe handling. Other sterilisation methods and cycles may also be suitable, however, individuals or hospitals are advised to validate whichever method they deem appropriate at their institution. Ethylene Oxide (EtO) sterilisation and cold sterilisation techniques are not recommended.

Corin TriFit TS Hip, showing the modular taper head and stem.

The information contained in this package insert was current on the date it was printed but the package insert may have been revised after that date. To obtain a current package insert, please contact Corin.

Corin TriFit TS Hip, showing the modular taper head and stem.

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Corin France SAS