


 **OMNIlife science, Inc.**
480 Paramount Dr.
Raynham, MA 02767 USA
www.omnils.com
(508)824-2444

Rx: Federal law (USA) restricts this device to sale or use by or on the order of a physician.

 Corin France SAS
157 Rue Lavoisier
38330 Montbonnot Saint-Martin
France



APEX Knee™ System and APEX Hip™ System X-Ray Templates

DESCRIPTION

APEX X-Ray Templates instruments are supplied by **OMNIlife science** specifically to assist in approximating the implant size of the Apex Knee™ or Apex Hip™ System implants prior to the surgical procedure.

INTENDED USE

The X-Ray Templates have been designed to estimate the implant size of the Apex Knee™ or Apex Hip™ Systems by placing these templates over the patient x-ray. A Surgical Technique can be obtained from **OMNIlife Science** or its representatives.

STORAGE

Store in a clean, dry and temperate place.

WARNINGS AND PRECAUTIONS

- **OMNIlife science** devices must only be used by surgeons who have been trained in the surgical technique and are familiar with the instruments or devices provided.
- Prior to using the x-ray templates and corresponding instrument systems, the surgeon should give careful consideration to all aspects of the surgical intervention as well as the limitations of the implant and instruments.

All rights reserved. Apex Knee™ System and Apex Hip™ System are trademarks of **OMNIlife science, Inc.** Additional information may be obtained from **OMNIlife science, Inc.**